



Instructions for Catalog # 1150

Air and Emissions Particulate Matter on Filter Paper

Revision 090119

Description:

- This standard is packaged in a 50 mm polystyrene Petri dish containing a single 47 mm glass fiber filter.
- This standard is designed for use with CARB Method 5 and EPA Methods 5, 5A, 5B, 5D, and 5F.
- This standard is not preserved.
- The standard can be stored at room temperature.
- Tare weights for the filter are provided for use in calculating your results.
- This product is intended to be used as a quality control check of the entire analytical process for the analytes/matrix included in the standard.
- This sample must be consumed in its entirety. Using a smaller sample size may invalidate the assigned value and/or uncertainty shown on the certificate of analysis.
- The certified values apply to the sample after following the stated instructions.

Helpful Hints:

- This standard should be analyzed following the procedure specified in CARB Method 5 or EPA Methods 5, 5A, 5B, 5D, and 5F.
- This standard may be desiccated and weighed or transferred to a tared weighing dish.

Instructions:

1. This standard is ready for analysis as received.
2. Analyze the standard following the procedure specified in CARB Method 5 or EPA Methods 5, 5A, 5B, 5D, and 5F.

Safety:

ERA products may be hazardous and are intended for use by professional laboratory personnel trained in the competent handling of such materials. Responsibility for the safe use of these products rests entirely with the buyer and/or user. Safety Data Sheets (SDS) for all ERA products are available through our website www.eraqc.com.